

Anesoft

Anesthesia & Critical Care Software
1 Montpellier Laguna Niguel, CA 92677 USA

Dear Anesoft Customer,

Thank you for purchasing Anesoft **Critical Care Simulator 3**. We hope you enjoy the simulated cases. The program will operate in any version of Windows.

Follow the installation instructions below to install Critical Care Simulator 3:

1. If you are installing from a CD, start My Computer, browse to the file setup_crit3.msi on the CD-ROM, and copy it to your desktop.
2. If you are installing from the download, SAVE the download file setup_crit3.msi to your desktop.
3. Close all applications and double-click setup_crit3.msi to start the installation.
4. Follow the installation instructions in the setup program except:

For Windows 7 and Windows 8, install the program to

c:\Anesoft\Critical Care Simulator 3
instead of c:\Program Files(x86)\Anesoft\Critical Care Simulator 3

5. If you are in a country that uses a comma as the decimal delimiter rather than a period, change your Language/Region settings in control panel since Critical Care Simulator 3 requires a period as the decimal delimiter.
6. The Critical Care Simulator 3 installation program will create a shortcut in your Programs list and on your desktop.

Note: The easiest way to learn to use the Critical Care Simulator 3 program is to follow the step-by-step instructions for case 1. Select Help > Instructions after you get the case started.

This program was planned and produced in accordance with ACCME essentials. You can earn AMA PRA Category 2 credits with Critical Care Simulator 3. Each physician should claim only those hours that he/she actually spent in the educational activity. Typically up to 3 credits are claimed for this activity (one half credit per case). Please email a request to jill@anesoft.com if you would like us to send you a CME certificate.

Please email us at jill@anesoft.com if you have any questions, problems or comments.

Best regards,
Jill, Anesoft Customer Service
jill@anesoft.com